

LAB TRAINING EUPHRESKO PHYTFIRE.

LISBON, 3-4 OCTOBER 2012.

Samples:

Set of samples for each of the following three bacterial pathogens:

E. amylovora (A), *E. pyrifoliae* (Py), *E. piriflorinigrans* (Pfn)

Set of samples	Code used on tubes	
S Bacterial suspensions of pure cultures in water	S-A S-Py S-Pfn	Bacterial concentration (cfu/ml): 10^6 , 10^5 , 10^4 , 10^3 , 10^2 , 10
P Plant material (pear leaves and shoots) extract spiked with bacterial suspensions	P-A P-Py P-Pfn	Bacterial concentration (cfu/ml) in plant material: 10^6 , 10^5 , 10^4 , 10^3 , 10^2 , 10
T DNA extractions of plant material (leaves and shoots) extract spiked with bacterial suspensions, according to Taylor <i>et al.</i> (2001)	T-A T-Py T-Pfn	DNA extractions of samples spiked at bacterial concentration (cfu/ml): 10^6 , 10^5 , 10^4 , 10^3 , 10^2 , 10

Controls:

Positive and negative controls could be:

Negative controls	Positive controls
1) Water used for the suspensions 2) Water used for PCR mix 3) Healthy plant material extract 4) Taylor extraction buffer	1) 10^6 cfu/ml bacterial suspension 2) Plant material spiked with 10^6 cfu/ml bacterial suspension 3) DNA extract from target organism

REAL TIME PCR for *E. amylovora* (Pirc et al., 2009)

Date _____

Components	x 1	
Ultrapure water	2,5 µl	
2X PCR Master Mix	12,5 µl	
Primer Ams 116F (10µM)	2,25 µl	
Primer Ams 189R (10µM)	2,25 µl	
Probe Ams 141T (10µM)	0,5 µl	
Aliquot to	20 µl	
Sample DNA	5 µl	

Samples:

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

PCR Conditions:

95° C 10 min

95° C 15 sec
60° C 1 min } 40 cycles

REAL TIME PCR for *E. amylovora* (Gottsberger *et al.*, 2010)

Date _____

Components	x 1	
Ultrapure water	7,95 µl	
2X PCR Master Mix	10 µl	
Primer hpEaF (10µM)	0,5 µl	
Primer hpEaR (10µM)	0,5 µl	
Probe hpEaP (10µM)	0,05 µl	
Aliquot to	19 µl	
Sample DNA	1 µl	

Samples:

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

PCR Conditions:

95° C 10 min

95° C 15 sec
60° C 1 min } 40 cycles

CONVENTIONAL PCR for *E. amylovora* (Taylor et al., 2001)

Date _____

Components	x 1	
Ultrapure water	14,3 µl /14,1 µl	
10X Thermopol buffer	2,5 µl	
MgCl ₂ (50 mM)	0,75 µl	
dNTPs (10 mM)	0,25 µl	
Primer G1-F (10µM)	1 µl	
Primer G2-R (10µM)	1 µl	
Taq polymerase (5U/µl)	0,2 µl /0,4 µl	
Aliquot to	20 µl	
Sample DNA	5 µl	

Samples:

1		11		21		31	
2		12		22		32	
3		13		23		33	
4		14		24		34	
5		15		25		35	
6		16		26		36	
7		17		27		37	
8		18		28		38	
9		19		29		39	
10		20		30		40	

PCR Conditions:

95° C 3 min

94° C	30 sec	} 40 cycles
60° C	30 sec	
72° C	1 min	

72° C 5 min

15 °C ∞

Expected product size: 187 bp Agarose gel (%): 2 %

LAMP PROTOCOL for *E. amylovora* (amsL gene)

Date _____

Components	X1	
10x Thermopol buffer	5 µl	
dNTPs (10 mM)	5 µl	
MgSO ₄ (100mM)	2 µl	
BSA (10 mg/ml)	2 µl	
ALB FIP (100 µM)	1,2 µl	
ALB BIP (100 µM)	1,2 µl	
ALB F (10 µM)	1 µl	
ALB B (10 µM)	1 µl	
Bst DNA polymerase (8U/(100 µl)	2 µl	
Ultrapure water	24,6 µl	
Sample DNA	5 µl	

Samples:

Samples	Results

LAMP Conditions:

65° C 60 min

CONVENTIONAL PCR for *E. pyrifoliae* (Wensing *et al.*, 2011)

Date _____

Components	x 1	
Ultrapure water	36,3 µl /36,1 µl	
10X Thermopol buffer	5 µl	
MgCl ₂ (50 mM)	1,5 µl	
dNTPs (10 mM)	1 µl	
Primer EPPSGL1646 (10µM)	0,5 µl	
Primer EPPSGL2698c (10µM)	0,5 µl	
Taq polymerase (5U/µl)	0,2 µl /0,4 µl	
Aliquot to	45 µl	
Sample DNA	5 µl	

Samples:

1		11		21		31	
2		12		22		32	
3		13		23		33	
4		14		24		34	
5		15		25		35	
6		16		26		36	
7		17		27		37	
8		18		28		38	
9		19		29		39	
10		20		30		40	

PCR Conditions:

95° C 3 min

94° C 30 sec
 57° C 30 sec
 72° C 1 min

} 40 cycles

72° C 5 min

15° C ∞

Expected product size: 1052 bp

Agarose gel (%): 2 %

CONVENTIONAL PCR for *E. pyrifoliae* (Kim et al., 2001)

Date _____

Components	x 1	
Ultrapure water	35,3 µl / 35,1 µl	
10X Thermopol buffer	5 µl	
MgCl ₂ (50 mM)	1,5 µl	
dNTPs (10 mM)	1 µl	
Primer CPS1 (10µM)	1 µl	
Primer CPS2 (10µM)	1 µl	
Taq polymerase (5U/µl)	0,2 µl / 0,4 µl	
Aliquot to	45 µl	
Sample DNA	5 µl	

Samples:

1		11		21		31	
2		12		22		32	
3		13		23		33	
4		14		24		34	
5		15		25		35	
6		16		26		36	
7		17		27		37	
8		18		28		38	
9		19		29		39	
10		20		30		40	

PCR Conditions:

95° C 3 min

94° C 30 sec
 57° C 30 sec
 72° C 1 min

} 40 cycles

72° C 5 min
 15° C ∞

Expected product size: 1200 bp
 Agarose gel (%): 2 %

REAL TIME PCR for *E. pyrifoliae* (based on Wensing *et al.*, 2011)

Date _____

Components	x 1	
Ultrapure water	7,25 µl	
2X PCR Master Mix	12,5 µl	
Primer EPPSGS1089Q (10µM)	1 µl	
Primer EPPSGS1228Qc (10µM)	1 µl	
Probe EPPSGS (10µM)	0,25 µl	
Aliquot to	22 µl	
Sample DNA	3 µl	

Samples:

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

PCR Conditions:

95° C 10 min

95° C 15 sec
60° C 1 min } 45 cycles

CONVENTIONAL PCR for *E. piriflorinigrans* (Barbé et al., 2012)

Date _____

Components	x 1	
Ultrapure water	34,4 µl /34,2 µl	
10X Thermopol buffer	5 µl	
MgCl ₂ (50 mM)	3 µl	
dNTPs (10 mM)	0,4 µl	
Primer P37F (10µM)	1 µl	
Primer P37R (10µM)	1 µl	
Taq polymerase (5U/µl)	0,2 µl /0,4 µl	
Aliquot to	45 µl	
Sample DNA	5 µl	

Samples:

1		11		21		31	
2		12		22		32	
3		13		23		33	
4		14		24		34	
5		15		25		35	
6		16		26		36	
7		17		27		37	
8		18		28		38	
9		19		29		39	
10		20		30		40	

PCR Conditions:

94° C 3 min

94° C	30 sec	} 40 cycles
56° C	30 sec	
72° C	45 sec	

72° C 10 min

8 ° C ∞

Expected product size: 255 bp Agarose gel (%): 2 %
--

REAL TIME PCR for *E. piriflorinigrans* (Barbé et al., 2012)

Date _____

Components	x 1	
Ultrapure water	1,2 µl /3,12 µl	
2X PCR Master Mix	6 µl /12,5 µl	
Primer p37F (10µM)	0,8 µl /1,95 µl	
Primer p37R (10µM)	0,8 µl /1,95 µl	
Probe p37 (10µM)	0,2 µl /0,48 µl	
Aliquot to	9 µl /20 µl	
Sample DNA	3 µl/5 µl	

Samples:

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

PCR Conditions:

95° C 10 min

95° C 15 sec
70° C 1 min } 40 cycles